



Herolab

Manufacturer of Life Science
Instruments since 1979



RobotCen

Robotic Centrifuge

For automation in the laboratory or in the production

RobotCen

- For the automated laboratory
- For the automated production



Catalog No.	Item		Voltage	Phases	Hz
21 13 050	RobotCen		200-240 V	1	50/60 Hz
21 13 650	RobotCen	gekühlt	220-240 V	1	50 Hz

- Suitable for Swing out and Fixed angle rotors
- Individual configuration according to your needs
- Computer control of the centrifuge via serial interface
- Exact positioning of the rotor guaranties friction-free interaction with the robotic arm
- On request, lid with two hatches or alternatively with a big hatch for Dual Arm Robots
- With or without refrigeration according to your needs
- Suitable for applications with bottles, tubes or microplates
- Maximum capacity: 4 x 500 ml and up to 4 x 1.000 ml in bigger models



Herolab

All products built by Herolab have the CE sign. Herolab has its own Quality Management System which is according to ISO 9001.



WEEE-Reg.-Nr.
DE66734561

Herolab reserves the right to change technical specifications without prior notice.

2021-01

Herolab GmbH Laborgeräte
Ludwig-Wagner-Str. 12
69168 Wiesloch
GERMANY
Tel.: +49 (0) 6222 – 5802-0
Fax: +49 (0) 6222 – 5802-34
E-Mail: info@herolab.de
Internet: www.herolab.de