

Centrifuge Tubes with Sealing Caps – Advanced Oak Ridge-Tubes



Advantages of the AOR-Tubes:

- The tubes are anti-twist and gas-tight through the unique combination of screw cap with conical sealing plug.
- The tubes can be used for High Speed requirements.
- The tubes are reusable.
- The dimensions correspond to the standard. So they fit into all rotors and adapters of well-known centrifuge manufacturers.

Cat. No. Tubes with Caps	Nominal Volume	Filling Volume	Dimensions (mm) Ø x Length	Material	Tested up to ... rpm	Tested up to ... x g	Pieces per Pack
25 24 01	10 ml	9 ml	16 x 80	PC	24.000	57.960	10
25 24 11	10 ml	8 ml	16 x 80	PPCO	24.000	57.960	10
25 24 41	10 ml	9 ml	16 x 80	PE	18.000	32.600	10
25 24 51	10 ml	9 ml	16 x 80	PF	20.000	40.250	2
25 29 51	30 ml	27 ml	25 x 92	PC	24.000	64.397	10
25 29 61	30 ml	25 ml	25 x 92	PPCO	24.000	64.397	10
25 29 65	30 ml	24 ml	25 x 92	PE	10.000	11.180	10
25 29 71	30 ml	26 ml	25 x 92	PF	24.000	64.397	2
25 32 01	50 ml	40 ml	29 x 103	PC	24.000	69.549	10
25 32 11	50 ml	38 ml	29 x 103	PPCO	24.000	69.549	10
25 32 15	50 ml	38 ml	29 x 103	PE	12.000	17.387	10
25 32 21	50 ml	42 ml	29 x 103	PF	18.000	39.121	2
25 32 41	80 ml	73 ml	38 x 105	PC	20.000	50.086	10
25 32 45	80 ml	76 ml	38 x 105	PPCO	20.000	50.086	10
25 32 47	80 ml	69 ml	38 x 105	PE	12.000	18.031	10
25 32 49	80 ml	73 ml	38 x 105	PF	15.000	28.174	2

Materials:

- PC Polycarbonate
- PPCO Polypropylene-Copolymer
- PE Polyethylene
- PF Polyflor